

Proteomic analysis of therapeutic biomarkers to guide treatment of patients with bone metastases

Fabiola Cecchi¹, Sarit Schwartz, Shankar Sellappan¹, Marlene Darfler¹, Kerry Scott¹, Sheeno Thyparambil¹, Michael Emmert-Buck², Yeoun Jin Kim¹, Todd Hembrough¹

¹NantOmics, Research and Development, Rockville, Maryland, USA, ²Cancer Center of Sarasota-Manatee, Sarasota, Florida, USA

BACKGROUND

- Personalized medicine has revolutionized cancer treatment, but treatment-related biomarkers can be challenging to assess in biopsies of metastatic bone lesions or bone cancer.
- Bone-containing biopsy tissues are routinely softened with rapid, strong-acid decalcification agents prior to sectioning; this can destroy cellular morphology, reduce protein immunoreactivity, and degrade the quality of DNA/RNA.
- We assessed the effect of acid decalcification on quantities of protein biomarkers expressed by archival tumor specimens from patients with various cancers.
- We quantified clinically relevant protein targets in decalcified bone biopsies from patients with various cancers using selected reaction monitoring (SRM).
- KRAS oncogenic mutation is an indicator that a tumor will not respond to EGFR-targeted therapy. We attempted to detect KRAS mutation at the protein level using data-independent acquisition (DIA).

METHODS




Figure 1. Liquid Tissue™ preparation followed by LC-SRM was used to quantify 27 biomarkers. LC-DIA was used to detect KRAS mutation.

RESULTS

Proteomics unaffected by decalcification conditions




Figure 2. (A) Prior to LM, tissue samples from 3 tumors were processed with Decal-Stat™ solution for 0, 1, 3, 12, and 24 hours (B) Total protein remained fairly stable with increasing decalcification time, in contrast to DNA. (C) Tumor expression of HER2 protein showed no relationship with decalcification time.

A clinically actionable oncogenic mutation was detected




Figure 3. Detection of mutant KRAS (G12D). DIA analysis surveying m/z [470-560] with 2Th window detected both wild type and mutant peptides in a bone biopsy.

SRM detected 18 of 27 treatment-related biomarkers tested in bony biopsies

Table 1. Proteins quantitated in biopsies of bone metastases or bone cancer (N=23)

Protein	Cutoff	Agent	Breast cancer			Lung cancer			Gastric			Sarcomas			Other		
			1	2	3	1	2	3	1	2	3	1	2	3	1	2	3
hENT1	100	gemcitabine	█	█	█	█	█	█	█	█	█	█	█	█	█	█	
TOP1	1340	irinotecan	█	█	█	█	█	█	█	█	█	█	█	█	█	█	
EGFR	1500	targeted	█	█	█	█	█	█	█	█	█	█	█	█	█	█	
TOP2A	1570	anthracycline	█	█	█	█	█	█	█	█	█	█	█	█	█	█	
HER2	750	targeted	█	█	█	█	█	█	█	█	█	█	█	█	█	█	
RRM1	700	gemcitabine	█	█	█	█	█	█	█	█	█	█	█	█	█	█	
ERCC1	75	platinum	█	█	█	█	█	█	█	█	█	█	█	█	█	█	
TUBB3	850	taxane	█	█	█	█	█	█	█	█	█	█	█	█	█	█	
KRAS	1650	targeted	█	█	█	█	█	█	█	█	█	█	█	█	█	█	
HER3	175	targeted	█	█	█	█	█	█	█	█	█	█	█	█	█	█	
MET	400	targeted	█	█	█	█	█	█	█	█	█	█	█	█	█	█	
MGMT	200	temozolomide	█	█	█	█	█	█	█	█	█	█	█	█	█	█	
AR	100	targeted	█	█	█	█	█	█	█	█	█	█	█	█	█	█	
AXL	100	targeted	█	█	█	█	█	█	█	█	█	█	█	█	█	█	
MSLN	500	targeted	█	█	█	█	█	█	█	█	█	█	█	█	█	█	
FGFR1-4	200	targeted	█	█	█	█	█	█	█	█	█	█	█	█	█	█	
FRalpha	1300	pemetrexed	█	█	█	█	█	█	█	█	█	█	█	█	█	█	
IGF1R	150	targeted	█	█	█	█	█	█	█	█	█	█	█	█	█	█	

█ Likely respond (> cutoff) █ Non actionable (< cutoff)
 █ Likely resistance (> cutoff) █ Not detected (< LOQ)

The 23 bone biopsy samples expressed 18 of the 30 protein targets tested. Half of samples expressed one or more protein markers of response to targeted therapy, and the vast majority expressed markers of response or resistance to conventional chemotherapies.

CONCLUSIONS

- A commonly used decalcifying solution had no discernable effects on mass spectrometric quantification of biomarker proteins in archived tumor samples.
- In decalcified bone biopsy specimens from cancer patients, an SRM assay quantitated 18 therapeutically relevant protein biomarkers that may inform selection of personalized cancer treatments.
- Mutant KRAS was detected at the protein level by DIA in a bone biopsy.
- Targeted proteomics can be used to rescue decalcified bone samples that are not evaluable by genomics due to insufficient quality or quantity of sequenceable DNA.